simeox
by PhysioAssist

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

of patients find a higher subjective effectiveness of therapy performance with Simeox vs CPT

of patients wish to pursue bronchial drainage with Simeox

This clinical study documented the risks and benefits of short-term use of Simeox in patients with cystic fibrosis. Simeox had a good safety profile and was feasible. The addition of the device to manual CPT has the potential to facilitate increased mucus clearance from the first drainage session.

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ULM STUDY

Simeox improves diaphragm mobility

Simeox therapy promotes the expectoration of highly viscous mucus from inflammatory areas

STUDY DESIGN

21 adults pwCF clinically stable; prospective, monocentric observational study
Comparison before/after 2 sessions of Simeox

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
Short term effects of a novel bronchial drainage device: a pilot cohort study in subjects with cystic fibrosis.

Health Sci Rep 2022;5:e812

doi.org/10.1002/hsr.2.812

Results


DIAPHRAGM EXCURSION


Maximum diaphragm excursion

SPUTUM VISCOSITY

Correlated with sputum DNA concentration

Sputum viscosity (FRAP)

Pearson correlation coefficient

3

●●● MAIN CLINICAL EVIDENCE FOR CYSTIC FIBROSIS

Patients at home

MOTHER AND CHILD INSTITUTE HOME CARE STUDY


MEF75 z-score (lung obstruction) improved significantly in Simeox group compared to Control group

Physical functioning score of CFQ-R improved significantly only in Simeox group


Results

LUNG FUNCTION


MEF75 z-score evolution from baseline at 1 month

QUALITY OF LIFE

Physical functioning score of CFQ-R evolution from baseline at 1 month

DISTAL AIRWAYS OBSTRUCTION

LC12.5 z-score evolution from baseline at 1 month

PATIENTS SURVEY

Recommend Simeox to CF patients

Prefer drainage with Simeox

Independent drainage

No fatigue

Relaxed exhalation

No drainage discomfort

● Definitively not ● Probably not ● Probably yes ● Definitively yes

Simeox improves safe mucus drainage of the airways in children with clinically stable CF under optimal therapy. Simeox better manages CF in children by decreasing airway obstruction in proximal airways, maintaining lung-clearance homogeneity, and increasing the CFQ-R physical functioning, all with safety and a high level of patient satisfaction.

LC12.5 z-score (ventilation heterogeneity) was preserved in Simeox group only after 1 month of therapy.

STUDY DESIGN

Comparison of safety and performance between Simeox plus optimal standard of care (SoC) and optimal SoC alone in 40 children with CF clinically stable (cross-over RCT, 2x1 month of therapy)

📍 Warsaw, Poland

REFERENCE

Sands, D.; Walicka-Serzysko, K.; Milczewska, J.; Postek, M.; Jeneralska, N.; Cichocka, A.; Siedlecka, E.; Borawska-Kowalczyk, U.; Morin, L. *Efficacy of the Simeox® Airway-Clearance Technology in the Homecare Treatment of Children with Clinically Stable Cystic Fibrosis: A Randomized Controlled Trial.*

Children 2023, 10, 204.

doi.org/10.3390/children10020204

On-going Multicenter RCTs

 IN CYSTIC FIBROSIS

MUCOSIM

Proof of concept (NCT04010253)

Based on positive outcomes of pilot studies in CF, this French multicenter, randomized controlled, non-inferiority trial was designed with the French Pulmonology Society to confirm in controlled environment the benefits and tolerance of one airway clearance session with Simeox in 35 clinically stable adult CF outpatients compared to autogenic drainage.

Primary endpoint is total airway resistance (R5hz) with IOS (Tremoflo®).

Recruitment is ongoing.

HEMOCARE CYSTIC FIBROSIS

(NCT04096664)

French multicenter, parallel group, superiority, randomized controlled trial designed to provide clinical evidence in the home setting on the efficacy and safety of Simeox in 160 clinically stable pwCF (>12y). The primary endpoints are in hierarchical order the Respiratory score of CFQ-R and FEV1 decline after 3 months of treatment in comparison with standard of care including usual CPT.

Remote setup and training of the device were successfully implemented since Covid-19 period.

Recruitment is ongoing.

Glossary

- **ACT:** Airway Clearance Techniques
- **CPT:** Chest Physiotherapy
- **CF:** Cystic Fibrosis
- **CFQ-R:** Cystic Fibrosis Questionnaire-Revised
- **CFTR:** Cystic Fibrosis Transmembrane conductance Regulator
- **DNA:** Deoxyribonucleic Acid
- **ENaC:** Epithelial sodium channel
- **FEV1:** Forced Expiratory Volume in 1 second
- **FRAP:** Fluorescence recovery after photobleaching
- **IOS:** Impulse Oscillometry System
- **LCI:** Lung Clearance Index
- **MEF25:** Maximal expiratory flow rate at 25% of vital capacity
- **MEF75:** Maximal expiratory flow rate at 75% of vital capacity
- **NS:** Not Significant
- **RCT:** Randomized Control Trial
- **PEx:** Pulmonary Exacerbation
- **pwCF:** Patients with Cystic Fibrosis
- **SoC:** Standard of Care

References

- (1) Flume PA, Robinson KA, O'Sullivan BP, Finder JD, Vender RL, Willey-Courand DB, White TB, Marshall BC, and the Clinical Practice Guidelines for Pulmonary Therapies Committee. Cystic Fibrosis Pulmonary Guidelines: Airway Clearance Therapies. *Respir Care* 2009;54(4): 522–537.
- (2) A.R. Smyth et al. European Cystic Fibrosis Society Standards of Care: Best Practice guideline. *Journal of Cystic Fibrosis* 13 (2014) S23–S42
doi.org/10.1016/j.jcf.2014.03.010.
- (3) Association of Chartered Physiotherapists in Cystic Fibrosis. *Standards of care and good clinical practice for the physiotherapy management of cystic fibrosis*. London: UK Cystic Fibrosis Trust; 2017.

PhysioAssist

31 Parc du Golf - CS 90519
13593 Aix-en-provence Cedex 3
FRANCE

tel +33 (0)4 67 03 13 92

contact@physio-assist.com

Simeox is a class IIa medical device and is a regulated healthcare product that carries, under the regulations, the CE mark. Read the instructions in the user manual carefully before using the device.

CMB_CF_EN01_0623